

# DK50 DS

## SMART

### SMALL AND COMPACT MEDICAL COMPRESSORS

#### ALARMS

- temperature ▣
- pressure ▣
- power ▣

#### INDICATIONS

- outlet pressure ▣
- working hours ▣
- drying ▣



Technical specification	DK50 DS	
Nominal voltage and frequency (*)	230/50, 230/60 115/50, 115/60 100/50..60	V/Hz
Compressor efficiency	40 (50***)	l.min <sup>-1</sup> / 3 bar
Peak flow at pressure 3 bar	200	l.min <sup>-1</sup>
Air tank capacity	2	lit.
Filtration of air to	5	µm
Dew point depression at 3 bar/ 40 lit.min <sup>-1</sup> , 20°C	5°C under the ambient temperature	
Dew point depression at 3 bar/ 40 lit.min <sup>-1</sup> , 20°C; with membrane dryer (**)	10°C under the ambient temperature	
Outlet connection	DISS 1160-A (3/4"-16 UNF)	
Mode of the operation	Continuous - S1	
Noise level	≤ 49	dB [A]
Dimensions (w x d x h)	445 x 355 x 440	mm
Weight	30	kg

\* please state voltage and frequency in the order

\*\* compressors can be completed with membrane dryer under the special order

\*\*\* under the special order